

## ***FLOWERS* Protocol**

### **Day 1**

- ☐ Make Stigma-Based Medium (adapted)
- ☐ Make starters of jdc\_isolates and *E. amylovora* in LB; incubate at 30°C

### **Day 2**

- ☐ To make standardized starters, transfer to Stigma-Based Medium:
    - ☐ jdc\_isolates starters (6µL to 600µL of SBM per well in 96 deep well plate)
    - ☐ *E. amylovora* starter (20µL to 2mL of SBM in snap cap Falcon tube)
- ...and incubate at 30°C

### **Day 3**

- ☐ Wait for standardized starters to grow

### **Day 4**

- ☐ Measure OD of jdc\_isolates and *E. amylovora* standardized starters
- ☐ Assemble communities using standardized starters of jdc\_isolates:
  - ☐ plate\_singles\_1 = jdc\_isolates
  - ☐ plate\_singles\_2 = jdc\_isolates.rotate\_left()
  - ☐ plate\_pairs\_1 = plate\_singles\_1 + plate\_singles\_2
  - ☐ plate\_pairs\_2 = plate\_pairs\_1.shuffle\_rows()
  - ☐ plate\_quad\_1 = plate\_pairs\_1 + plate\_pairs\_2
  - ☐ plate\_quad\_2 = plate\_quad\_1.shuffle\_rows()
  - ☐ plate\_quad\_3 = plate\_quad\_2.rotate\_left()
  - ☐ plate\_communities = plate\_quad\_1 + plate\_quad\_3
- ☐ Inoculate communities and standardized starter of *E. amylovora* into fresh SBM (2µL of each into 200µL of SBM per well in 96 deep well plate)
- ☐ In a different control plate, inoculate communities into fresh SBM (2µL into 200µL of SBM per well in 96 deep well plate)
- ☐ Incubate at 30°C with transpirable filter

### **Day 5**

- ☐ Wait for communities to grow

### **Day 6**

- ☐ Wait for communities to grow

**Day 7**

- ☐ Transfer 100 $\mu$ L per well to a standard 96 well plate
- ☐ Measure OD
- ☐ Measure GFP fluorescence
- ☐ Use FACS to quantify abundance of *E. amylovora* in each community
  - ☐ If needed: dilute 30 $\mu$ L of the communities into 270 $\mu$ L of PBS per well in 96 deep well plate (1:10 dilution for a final 10-100 million cells per mL as required by FACS specifications)