

***FLOWERS* Protocol**

Day 1

- ☐ Make Stigma-Based Medium (adapted)
- ☐ Make starters of jdc_isolates and *E. amylovora* in SBM; incubate at 30°C

Day 2

- ☐ Wait for starters to grow

Day 3

- ☐ To make standardized starters, transfer to Stigma-Based Medium:
 - ☐ jdc_isolates starters (6µL to 600µL of SBM per well in 96 deep well plate)
 - ☐ *E. amylovora* starter (20µL to 2mL of SBM in snap cap Falcon tube)

Day 4

- ☐ Wait for standardized starters to grow

Day 5

- ☐ Measure OD of jdc_isolates and *E. amylovora*
- ☐ Assemble communities using standardized starters of jdc_isolates:
 - ☐ plate_singles_1 = jdc_isolates
 - ☐ plate_singles_2 = jdc_isolates.rotate_left()
 - ☐ plate_pairs_1 = plate_singles_1 + plate_singles_2
 - ☐ plate_pairs_2 = plate_pairs_1.shuffle_rows()
 - ☐ plate_quad_1 = plate_pairs_1 + plate_pairs_2
 - ☐ plate_quad_2 = plate_quad_1.shuffle_rows()
 - ☐ plate_quad_3 = plate_quad_2.rotate_left()
 - ☐ plate_communities = plate_quad_1 + plate_quad_3
- ☐ Inoculate communities and standardized starter of *E. amylovora* into fresh SBM (2µL of each into 200µL of SBM per well in 96 deep well plate)
- ☐ Incubate at 30°C with transpirable filter

Day 6

- ☐ Wait for communities to grow

Day 7

- ☐ Wait for communities to grow

Day 8

- ☐ Transfer 100 μ L per well to a standard 96 well plate
- ☐ Measure OD
- ☐ Measure GFP fluorescence
- ☐ Use FACS to quantify abundance of *E. amylovora* in each community
 - ☐ If needed: dilute 30 μ L of the communities into 270 μ L of PBS per well in 96 deep well plate (1:10 dilution for a final 10-100 million cells per mL as required by FACS specifications)